bs-9903R

[Primary Antibody]

SCUBE1 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GeneID: 80274 SWISS: Q8IWY4

Target: SCUBE1

Immunogen: KLH conjugated synthetic peptide derived from human SCUBE1:

451-550/988.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: Could function as an adhesive molecule and its matrix bound and

soluble fragments may play a critical role in vascular biology.

Applications: WB (1:500-2000)

IHC-P (1:100-500) **IHC-F** (1:100-500) **IF** (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit, Pig, Cow, Dog,

Horse)

Predicted MW.: 105 kDa

Subcellular Location: Secreted ,Cell membrane