

**bs-0373R-PE-Cy7**

**[ Secondary Antibodies ]**

## **Rabbit Anti-Dog IgM, PE-Cy7 conjugated**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> IF (1:100-1000)
<b>Clonality:</b> Polyclonal		<b>Flow-Cyt</b> (1:100-1000)
<b>Target:</b> Rabbit Anti-Dog IgM		<b>ICC/IF</b> (1:100-1000)
<b>Purification:</b> affinity purified by Protein A		Excitation Spectrum: 488nm
<b>Concentration:</b> 2.0 mg/ml		Emission spectrum: 780nm
<b>Storage:</b> 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		<b>Reactivity:</b> Dog