

**bs-9082R**

**[ Primary Antibody ]**

## AKIRIN2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 55122</p> <p><b>Target:</b> AKIRIN2</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human AKIRIN2: 51-150/203.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> Required for the innate immune response. Downstream effector of the Toll-like receptor (TLR), TNF and IL-1 beta signaling pathways leading to the production of IL-6. Forms a complex with YWHAB that acts to repress transcription of DUSP1.</p>	<p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q53H80</p>	<p><b>Applications:</b> <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:50-200) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> Mouse (predicted: Human, Rat, Pig, Sheep, Cow, Chicken, Dog, Horse)</p> <p><b>Predicted MW.:</b> 22 kDa</p> <p><b>Subcellular Location:</b> Nucleus</p>
--	--	--