bs-9082R

[Primary Antibody]

Bioss ANTIBODIES

AKIRIN2 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 55122 **SWISS:** Q53H80

Target: AKIRIN2

Immunogen: KLH conjugated synthetic peptide derived from human AKIRIN2:

51-150/203.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: Required for the innate immune response. Downstream effector of

the Toll-like receptor (TLR), TNF and IL-1 beta signaling pathways leading to the production of IL-6. Forms a complex with YWHAB

that acts to repress transcription of DUSP1.

Applications: IHC-P (1:100-500)

400-901-9800

IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: Mouse (predicted: Human,

Rat, Pig, Sheep, Cow, Chicken, Dog, Horse)

Predicted MW.: 22 kDa

Subcellular Nucleus