

---

## RBM24 Rabbit pAb

Catalog Number: bs-21293R

Target Protein: RBM24

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse)

Predicted MW: 25 kDa

Entrez Gene: 221662

Swiss Prot: Q9BX46

Source: KLH conjugated synthetic peptide derived from human RBM24: 31-130/236.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: RBM24 (RNA Binding Motif Protein 24) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and mRNA 3-UTR binding. An important paralog of this gene is RBM38.