

bs-21293R

[**Primary Antibody**]

RBM24 Rabbit pAb

BioSS[™]
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 221662</p> <p>Target: RBM24</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human RBM24: 31-130/236.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: RBM24 (RNA Binding Motif Protein 24) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and mRNA 3-UTR binding. An important paralog of this gene is RBM38.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9BX46</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse)</p> <p>Predicted MW.: 25 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
---	--	--