
phospho-Tau (Thr231) Mouse mAb

Cat. No:	V5221
Name:	phospho-Tau (Thr231) Mouse mAb
Clone type:	Monoclonal
Clone No:	7F12
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr231 (bs-25325P)