

**bs-2146R**

**[ Primary Antibody ]**

## HBV pre-X protein Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|  |  |
|--|--|
| <p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> HBV pre-X protein</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from Hepatitis B virus X protein: 1-50/210.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br/>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b></p> | <p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Hepatitis B virus)</p> <p><b>Predicted MW.:</b> 23 kDa</p> <p><b>Subcellular Location:</b> Cell membrane ,Cytoplasm</p> |
|--|--|