

---

## ZNF430 Rabbit pAb

Catalog Number: bs-7680R

Target Protein: ZNF430

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 66 kDa

Entrez Gene: 80264

Swiss Prot: Q9H8G1

Source: KLH conjugated synthetic peptide derived from human ZNF430: 671-570/570.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.