

bs-4456R

[Primary Antibody]

ZG16A Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 653808</p> <p>Target: ZG16A</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZG16A: 31-120/167.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZG16 (Zymogen Granule Protein 16) is a Protein Coding gene. GO annotations related to this gene include carbohydrate binding. An important paralog of this gene is ZG16B.</p>	<p>Isotype: IgG</p> <p>SWISS: O60844</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Horse)</p> <p>Predicted MW.: 18 kDa</p> <p>Subcellular Location: Extracellular matrix</p>
---	--	---