

**bs-8749R**

**[ Primary Antibody ]**

## THSD7A Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|                                                                                                                                                                      |                      |                                                                                                                                                                                                                                                                                                                                             |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host:</b> Rabbit                                                                                                                                                  | <b>Isotype:</b> IgG  | <b>Applications:</b> <b>IHC-P</b> (1:100-500)<br><b>IHC-F</b> (1:100-500)<br><b>IF</b> (1:100-500)<br><b>ICC/IF</b> (1:100-500)<br><b>ELISA</b> (1:5000-10000)<br><br><b>Reactivity:</b> (predicted: Human, Mouse, Rat, Rabbit, Chicken, Dog, Horse)<br><br><b>Predicted MW.:</b> 180 kDa<br><br><b>Subcellular Location:</b> Cell membrane |
| <b>Clonality:</b> Polyclonal                                                                                                                                         |                      |                                                                                                                                                                                                                                                                                                                                             |
| <b>GeneID:</b> 221981                                                                                                                                                | <b>SWISS:</b> Q9UPZ6 |                                                                                                                                                                                                                                                                                                                                             |
| <b>Target:</b> THSD7A                                                                                                                                                |                      |                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen:</b> KLH conjugated synthetic peptide derived from human THSD7A: 1201-1300/1657. < Extracellular >                                                      |                      |                                                                                                                                                                                                                                                                                                                                             |
| <b>Purification:</b> affinity purified by Protein A                                                                                                                  |                      |                                                                                                                                                                                                                                                                                                                                             |
| <b>Concentration:</b> 1mg/ml                                                                                                                                         |                      |                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |                      |                                                                                                                                                                                                                                                                                                                                             |