

bs-12086R

[Primary Antibody]

JAKMIP1/GABABRBP Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 152789</p> <p>Target: JAKMIP1/GABABRBP</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human GABABRBP/JAKMIP1: 20-120/626.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Associates with microtubules and may play a role in the microtubule-dependent transport of the GABA-B receptor. May play a role in JAK1 signaling and regulate microtubule cytoskeleton rearrangements.</p>	<p>Isotype: IgG</p> <p>SWISS: Q96N16</p> <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 73 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	---