
JAKMIP1/GABABRBP Rabbit pAb

Catalog Number: bs-12086R

Target Protein: JAKMIP1/GABABRBP

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 73 kDa

Entrez Gene: 152789

Swiss Prot: Q96N16

Source: KLH conjugated synthetic peptide derived from human GABABRBP/JAKMIP1: 20-120/626.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Associates with microtubules and may play a role in the microtubule-dependent transport of the GABA-B receptor. May play a role in JAK1 signaling and regulate microtubule cytoskeleton rearrangements.