

bs-1690R

[Primary Antibody]

H2N2 Hemagglutinin Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: H2N2 Hemagglutinin

Immunogen: KLH conjugated synthetic peptide derived from H2N2 Hemagglutinin: 461-562/562.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: No data available.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Influenza A virus H2N2)

Predicted MW.: 62 kDa