

**bs-12598R**

**[ Primary Antibody ]**

## PDILT Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|  |                      |   |
|--|----------------------|---|
| <b>Host:</b> Rabbit  | <b>Isotype:</b> IgG  | <b>Applications:</b> <b>IHC-P</b> (1:100-500)<br><b>IHC-F</b> (1:100-500)<br><b>IF</b> (1:100-500)<br><b>ICC/IF</b> (1:100-500)<br><b>ELISA</b> (1:5000-10000)<br><br><b>Reactivity:</b> (predicted: Human, Mouse, Rat)<br><br><b>Predicted MW.:</b> 64 kDa<br><br><b>Subcellular Location:</b> Cytoplasm |
| <b>Clonality:</b> Polyclonal   |                      |   |
| <b>GeneID:</b> 204474  | <b>SWISS:</b> Q8N807 |   |
| <b>Target:</b> PDILT   |                      |   |
| <b>Immunogen:</b> KLH conjugated synthetic peptide derived from human PDILT: 201-300/584.  |                      |   |
| <b>Purification:</b> affinity purified by Protein A  |                      |   |
| <b>Concentration:</b> 1mg/ml   |                      |   |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |                      |   |