
SFMBT2 Rabbit pAb

Catalog Number: bs-21105R

Target Protein: SFMBT2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Cow, Horse)

Predicted MW: 101 kDa

Entrez Gene: 57713

Swiss Prot: Q5VUG0

Source: KLH conjugated synthetic peptide derived from human SFMBT2: 201-300/894.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: SFMBT2 (Scm-Like With Four Mbt Domains 2) is a Protein Coding gene. GO annotations related to this gene include sequence-specific DNA binding and histone binding. An important paralog of this gene is SCMH1.