

bs-21105R

[Primary Antibody]

SFMBT2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 57713</p> <p>Target: SFMBT2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SFMBT2: 201-300/894.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: SFMBT2 (Scm-Like With Four Mbt Domains 2) is a Protein Coding gene. GO annotations related to this gene include sequence-specific DNA binding and histone binding. An important paralog of this gene is SCMH1.</p>	<p>Isotype: IgG</p> <p>SWISS: Q5VUG0</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Cow, Horse)</p> <p>Predicted MW.: 101 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	---