

**bs-8530R**

**[ Primary Antibody ]**

## SDHC Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 6391</p> <p><b>Target:</b> SDHC</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human SDHC: 81-169/169.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> Membrane-anchoring subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone.</p>	<p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q99643</p>	<p><b>Applications:</b> <b>WB</b> (1:500-2000) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human)</p> <p><b>Predicted MW.:</b> 15 kDa</p> <p><b>Subcellular Location:</b> Cell membrane ,Cytoplasm</p>
--	--	--