

bs-37015R

[Primary Antibody]

Enfuvirtide Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Enfuvirtide

Purification: affinity purified by Protein A

Concentration: 1mg/1ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Enfuvirtide (INN) is an HIV fusion inhibitor, the first of a novel class of antiretroviral drugs used in combination therapy for the treatment of HIV-1 infection. It is marketed under the trade name Fuzeon (Roche).

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: Enfuvirtide)

**Predicted
MW.:** 4.49 kDa