

**bs-37015R**

**[ Primary Antibody ]**

## Enfuvirtide Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Enfuvirtide

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/1ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Enfuvirtide (INN) is an HIV fusion inhibitor, the first of a novel class of antiretroviral drugs used in combination therapy for the treatment of HIV-1 infection. It is marketed under the trade name Fuzeon (Roche).

**Applications:** ELISA (1:5000-10000)

**Reactivity:** Species independent

**Predicted MW.:** 4.49 kDa