

## phospho-Tau (Thr231) Mouse mAb

|                      |                                                                                 |
|----------------------|---------------------------------------------------------------------------------|
| Cat. No:             | V5222                                                                           |
| Name:                | phospho-Tau (Thr231) Mouse mAb                                                  |
| Storage:             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:          | Monoclonal                                                                      |
| Clone No:            | 4C10                                                                            |
| Host/Source:         | Mouse                                                                           |
| Subtype:             | IgG2b                                                                           |
| Reactivity:          | (predicted:&nbsp;Human)                                                         |
| Purity:              | >90% as determined by SDS-PAGE                                                  |
| Format:              | Liquid                                                                          |
| Concentration:       | 1mg/ml                                                                          |
| Buffer:              | 0.01M PBS (pH7.4).                                                              |
| Applications:        | ELISA                                                                           |
| Antigen:             | Recombinant Human p-Tau (Thr231) Protein ( <a href="#">bs-41576P</a> )          |
| Pairing Information: |                                                                                 |

| ELISA   |       | Detection |       |       |       |       |       |       |       |       |       |       |       |       |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
|         |       | V5203     | V5206 | V5214 | V5215 | V5216 | V5218 | V5219 | V5232 | V5233 | V5234 | V5236 | V5237 | W5201 | V5209 | V5210 | V5211 | V5226 |
| Capture | V5203 | -         | -     | +     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5204 | -         | -     | +     | -     | +     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5205 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5206 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5207 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5212 | +         | -     | -     | -     | -     | -     | +     | +     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5213 | +         | +     | +     | +     | +     | -     | -     | -     | -     | -     | -     | -     | -     | +     | +     | +     | -     |
|         | V5215 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     |
|         | V5216 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5217 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     |
|         | V5218 | -         | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     |
|         | V5220 | -         | -     | -     | -     | -     | -     | -     | -     | -     | +     | -     | -     | +     | -     | -     | -     | -     |
|         | V5231 | -         | -     | -     | -     | -     | -     | -     | -     | +     | +     | -     | -     | -     | -     | -     | -     | -     |
|         | V5233 | -         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | +++   |
|         | V5234 | -         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | +     |
|         | V5235 | -         | -     | -     | -     | -     | -     | +     | +     | -     | -     | +     | -     | -     | -     | -     | -     | -     |
|         | V5236 | -         | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     | -     | -     | -     | -     | +++   |
|         | W5201 | -         | -     | -     | -     | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     | -     | -     |
|         | V5225 | -         | +     | +     | -     | -     | -     | +     | -     | -     | -     | +     | -     | -     | -     | -     | -     | -     |
|         | V5227 | +         | -     | -     | -     | -     | -     | -     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     |
| V5222   | -     | -         | +     | -     | -     | -     | +     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     |       |
| V5223   | -     | -         | +     | -     | -     | -     | +     | -     | -     | -     | -     | +     | -     | -     | -     | -     | -     |       |