

---

## phospho-Tau (Thr231) Mouse mAb

Cat. No:	V5224
Name:	phospho-Tau (Thr231) Mouse mAb
Clone type:	Monoclonal
Clone No:	7B3
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr231 ( <a href="#">bs-25325P</a> )