

**bs-20822R**

**[ Primary Antibody ]**

## TMEM120A Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 215210</p> <p><b>Target:</b> TMEM120A</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from mouse TMEM120A: 251-343/343. &lt; Extracellular &gt;</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q8C1E7</p>	<p><b>Applications:</b> <b>WB</b> (1:500-2000) <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ICC/IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Mouse, Rat, Sheep, Cow, Dog)</p> <p><b>Predicted MW.:</b> 41 kDa</p> <p><b>Subcellular Location:</b> Cell membrane</p>
--	--	---