

bs-21121R

[Primary Antibody]

QPCTL Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 54814 Target: QPCTL Immunogen: KLH conjugated synthetic peptide derived from human QPCTL: 181-280/382. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG SWISS: Q9NXS2	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Horse) Predicted MW.: 43 kDa Subcellular Location: Cytoplasm
---	---	---