

bs-21123R

[Primary Antibody]

QRSL Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 55278</p> <p>Target: QRSL</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human QRSL: 401-500/528.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: QRSL belongs to the amidase family. Catalytic activity: ATP + L-glutamyl-tRNA(Gln) + L-glutamine = ADP + phosphate + L-glutamyl-tRNA(Gln) + L-glutamate. There are two named isoforms.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9H0R6</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog)</p> <p>Predicted MW.: 57 kDa</p> <p>Subcellular Location: Cytoplasm</p>
--	--	--