

bs-18709R

[Primary Antibody]

MBOAT2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Chicken, Dog, Horse) Predicted MW.: 59 kDa Subcellular Location: Cell membrane
Clonality: Polyclonal		
GeneID: 129642	SWISS: Q6ZWT7	
Target: MBOAT2		
Immunogen: KLH conjugated synthetic peptide derived from human MBOAT2: 81-180/520.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=7.1]** Wei Deng. et al. Role of ferroptosis mediated by abnormal membrane structure in DEHP-induced reproductive injury. FREE RADICAL BIO MED. 2025 Aug;235:150 IHC ;Mouse. 40306442
- **[IF=3.9]** Bo Yu. et al. TMAO induces pyroptosis of vascular endothelial cells and atherosclerosis in ApoE—/— mice via MBOAT2-mediated endoplasmic reticulum stress. BBA-MOL CELL BIOL L. 2024 Aug;;159559 IF ;Mouse. 39179098