

bs-11187R

[Primary Antibody]

Krt42 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 68239 Target: Krt42 Immunogen: KLH conjugated synthetic peptide derived from mouse Krt42: 301-400/452. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: Predicted to enable structural molecule activity. Predicted to be located in cytoplasm and intermediate filament. [provided by Alliance of Genome Resources, Nov 2021]	Isotype: IgG SWISS: Q6IFX2	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: Human (predicted: Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse) Predicted MW.: 50 kDa Subcellular Location: Cytoplasm
--	--	--