

**bs-19338R**

**[ Primary Antibody ]**

## NPTXR Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> ELISA (1:5000-10000)
<b>Clonality:</b> Polyclonal		<b>Reactivity:</b> (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Horse, Chimpanzee)
<b>GeneID:</b> 23467	<b>SWISS:</b> O95502	
<b>Target:</b> NPTXR		<b>Predicted MW.:</b> 53 kDa
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human NPTXR: 221-320/500. < Extracellular >		<b>Subcellular Location:</b> Cell membrane
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		