

bs-18139R

[Primary Antibody]

IGKV1-5 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 28299 Target: IGKV1-5 Immunogen: KLH conjugated synthetic peptide derived from human IGKV1-5: 23-100/117. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: immunoglobulin kappa variable 1-5	Isotype: IgG SWISS: P01602	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 10 kDa Subcellular Location: Cell membrane
--	--	--