

**bs-21161R**

**[ Primary Antibody ]**

## RASAL3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 64926</p> <p><b>Target:</b> RASAL3</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human RASAL3: 721-820/1011.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q86YV0</p>	<p><b>Applications:</b> <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ICC/IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human, Mouse, Rat, Sheep, Cow, Dog)</p> <p><b>Predicted MW.:</b> 112 kDa</p> <p><b>Subcellular Location:</b> Cytoplasm</p>
--	--	--