

bs-21162R

[Primary Antibody]

RASGEF1B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 153020</p> <p>Target: RASGEF1B</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human RASGEF1B: 1-100/473.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: RASGEF1B (RasGEF Domain Family Member 1B) is a Protein Coding gene. Diseases associated with RASGEF1B include Chromosome 4Q21 Deletion Syndrome. GO annotations related to this gene include guanyl-nucleotide exchange factor activity and Ras guanyl-nucleotide exchange factor activity. An important paralog of this gene is RASGEF1C.</p>	<p>Isotype: IgG</p> <p>SWISS: Q0VAM2</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Pig, Sheep, Cow)</p> <p>Predicted MW.: 55 kDa</p> <p>Subcellular Location: Cytoplasm</p>
---	--	--