

---

## RASGEF1B Rabbit pAb

Catalog Number: bs-21162R

Target Protein: RASGEF1B

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Pig, Sheep, Cow)

Predicted MW: 55 kDa

Entrez Gene: 153020

Swiss Prot: Q0VAM2

Source: KLH conjugated synthetic peptide derived from human RASGEF1B: 1-100/473.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** RASGEF1B (RasGEF Domain Family Member 1B) is a Protein Coding gene. Diseases associated with RASGEF1B include Chromosome 4Q21 Deletion Syndrome. GO annotations related to this gene include guanyl-nucleotide exchange factor activity and Ras guanyl-nucleotide exchange factor activity. An important paralog of this gene is RASGEF1C.