

bs-21166R

[Primary Antibody]

SGOL2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 151246</p> <p>Target: SGOL2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SGOL2: 781-880/1265.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: SGO2 (Shugoshin 2) is a Protein Coding gene. Diseases associated with SGO2 include Perrault Syndrome. Among its related pathways are Mitotic Metaphase and Anaphase and Signaling by GPCR.</p>	<p>Isotype: IgG</p> <p>SWISS: Q562F6</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 145 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	---