

**bsm-4537M**

**[ Primary Antibody ]**

## **Streptococcus suis Suilysin Mouse mAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

<b>Host:</b> Mouse	<b>Isotype:</b> IgG	<b>Applications:</b> ELISA (1:5000-10000)  <b>Reactivity:</b> (predicted: Streptococcus suis Suilysin)
<b>Clonality:</b> Monoclonal	<b>CloneNo.:</b> 2H3	
<b>Target:</b> Streptococcus suis Suilysin		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> Size : 50ul/100ul/200ul Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		