

bs-12608R

[Primary Antibody]

PDZRN4 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 29951</p> <p>Target: PDZRN4</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PDZRN4: 761-860/1036.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: PDZRN4 (PDZ domain-containing RING finger protein 4) contains two PDZ (DHR) domains and a RING-type zinc finger. Its function is unknown.</p>	<p>Isotype: IgG</p> <p>SWISS: Q6ZMN7</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 117 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
--	--	--