
PDZRN4 Rabbit pAb

Catalog Number: bs-12608R

Target Protein: PDZRN4

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 117 kDa

Entrez Gene: 29951

Swiss Prot: Q6ZMN7

Source: KLH conjugated synthetic peptide derived from human PDZRN4: 761-860/1036.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: PDZRN4 (PDZ domain-containing RING finger protein 4) contains two PDZ (DHR) domains and a RING-type zinc finger. Its function is unknown.