

bs-21190R

[Primary Antibody]

SLC12A8 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 84561</p> <p>Target: SLC12A8</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SLC12A8: 501-600/714.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene is thought to be a candidate for psoriasis susceptibility. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Sep 2010] Function : Cation/chloride cotransporter that may play a role in the control of keratinocyte proliferation. Subcellular Location : Membrane.</p>	<p>Isotype: IgG</p> <p>SWISS: A0AV02</p>	<p>Applications: WB (1:500-2000) ELISA (1:5000-10000)</p> <p>Reactivity: Human</p> <p>Predicted MW.: 78 kDa</p> <p>Subcellular Location: Cell membrane</p>
--	--	--