

---

## SLC12A8 Rabbit pAb

Catalog Number: bs-21190R

Target Protein: SLC12A8

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), ELISA (1:5000-10000)

Reactivity: Human

Predicted MW: 78 kDa

Subcellular Cell membrane

Locations:

Entrez Gene: 84561

Swiss Prot: A0AV02

Source: KLH conjugated synthetic peptide derived from human SLC12A8: 501-600/714.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** This gene is thought to be a candidate for psoriasis susceptibility. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Sep 2010] Function : Cation/chloride cotransporter that may play a role in the control of keratinocyte proliferation. Subcellular Location : Membrane.