

bs-17456R

[Primary Antibody]

HIV1 p24 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: HIV1 p24 Immunogen: KLH conjugated synthetic peptide derived from HIV1 p24: 251-350/497. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: HIV1) Predicted MW.: 26 kDa
--	---

— SELECTED CITATIONS —

- **[IF=6.1]** Songbai Tian. et al. Membrane-based microfluidic chip constructed and applied in magnetic microparticle anchoring-and-release digital immunoassay. TALANTA. 2026 Jan;297:128644 ;. 40743937
- **[IF=5.256]** Du MY et al. Self-Assembled Fluorescent Ce (III) Coordination Polymer as Ratiometric Probe for HIV Antigen Detection. Analytica Chimica Acta. 2019. ELISA ;Human. doi:10.1016/j.aca.2019.08.010