

bsm-41186M**[Primary Antibody]****NGAL/Lipocalin 2(5B3) Mouse mAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Monoclonal	CloneNo.: 5B3	Reactivity: (predicted: Human)
GeneID: 3934	SWISS: P80188	
Target: NGAL/Lipocalin 2(5B3)		
Purification: affinity purified by Protein G		Predicted MW.: 20 kDa
Concentration: 1mg/ml		Subcellular Location: Secreted
Storage: Size : 50ul/100ul/200ul 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Lipocalin 2 is a member of the lipocalin family which encompass more than 25 members (including proteins like retinol binding protein, a1 microglobulin, b lactoglobulin, apolipoprotein D and odorant binding protein). The members of the lipocalin family are characterized by their ability to bind small lipophilic substances in their hydrophobic core, and thereby serve as transporters of substances.		