## bs-21204R

## [ Primary Antibody ]

## Bioss

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **SLC1A7 Rabbit pAb**

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**GenelD:** 6512 **SWISS:** 000341

Target: SLC1A7

**Immunogen:** KLH conjugated synthetic peptide derived from human SLC1A7:

1-100/560. < Cytoplasmic >

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS,

pH1.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)
IF (1:100-500)
ICC/IF (1:100-500)
ELISA (1:5000-10000)

Reactivity: (predicted: Human, Sheep,

Cow, Dog, Horse)

Predicted MW.: 61 kDa

Subcellular Location: Cell membrane

\_\_\_\_\_\_