
Histone H2A.V Rabbit pAb

Catalog Number: bs-17425R

Target Protein: Histone H2A.V

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Cow, Zebrafish, Chicken, Dog, Horse, Xenopus laevis)

Predicted MW: 14 kDa

Subcellular Nucleus

Locations:

Entrez Gene: 94239

Swiss Prot: Q71UI9

Source: KLH conjugated synthetic peptide derived from human Histone H2A.V: 21-80/128.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

PRODUCT SPECIFIC PUBLICATIONS

[IF=1.664] Zhang PF et al. MicroRNA-139 suppresses hepatocellular carcinoma cell proliferation and migration by directly targeting Topoisomerase I. ONCOLOGY LETTERS 17: 1903-1913, 2019 WB ; Human . 10.3892/ol.2018.9746