

bs-17425R

[Primary Antibody]

Histone H2A.V Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Pig, Cow, Zebrafish, Chicken, Dog, Horse, Xenopus laevis) Predicted MW.: 14 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 94239	SWISS: Q71UI9	
Target: Histone H2A.V		
Immunogen: KLH conjugated synthetic peptide derived from human Histone H2A.V: 21-80/128.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=1.664]** Zhang PF et al. MicroRNA-139 suppresses hepatocellular carcinoma cell proliferation and migration by directly targeting Topoisomerase I. ONCOLOGY LETTERS 17: 1903-1913, 2019 WB ;Human. 10.3892/ol.2018.9746