

---

## ZNF576 Rabbit pAb

Catalog Number: bs-19525R

Target Protein: ZNF576

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Pig, Sheep, Cow, Horse)

Predicted MW: 19 kDa

Entrez Gene: 79177

Swiss Prot: Q9H609

Source: KLH conjugated synthetic peptide derived from human ZNF576: 1-100/170.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.