

bs-16482R

[Primary Antibody]

ZNF77 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 58492</p> <p>Target: ZNF77</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZNF77: 451-545/545.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZNF77 belongs to the krueppel C2H2 type zinc finger protein family. It contains 12 C2H2 type zinc fingers and 1 KRAB domain. ZNF77 may be involved in transcriptional regulation.</p>	<p>Isotype: IgG</p> <p>SWISS: Q15935</p> <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Horse)</p> <p>Predicted MW.: 62 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--