

**bs-0297G-BF594**

**[ Secondary Antibodies ]**

## **Goat Anti-Human IgG H&L, BF594 conjugated**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Goat

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Goat Anti-Human IgG H&L

**Purification:** affinity purified by Protein G, nonspecific adsorbed

**Concentration:** 2.0 mg/ml

**Storage:** 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.  
Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** IF (1:100-500)

**Flow-Cyt** (1:50-200)

**ICC/IF** (1:100-500)

**Reactivity:** Human