

bs-0297G-BF594

[Secondary Antibodies]

Goat Anti-Human IgG H&L, BF594 conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat

Isotype: IgG

Clonality: Polyclonal

Target: Goat Anti-Human IgG H&L

Purification: affinity purified by Protein G, nonspecific adsorbed

Concentration: 2.0 mg/ml

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.
Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IF (1:100-500)

Flow-Cyt (1:50-200)

ICC/IF (1:100-500)

Excitation Spectrum: 594nm

Emission spectrum: 621nm

Reactivity: Human