

bs-18721R

[Primary Antibody]

MCEMP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 21 kDa Subcellular Location: Cell membrane
Clonality: Polyclonal		
GeneID: 199675	SWISS: Q8IX19	
Target: MCEMP1		
Immunogen: KLH conjugated synthetic peptide derived from human MCEMP1: 101-187/187. < Extracellular >		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: This gene encodes a single-pass transmembrane protein. Based on its expression pattern, it is speculated to be involved in regulating mast cell differentiation or immune responses. [provided by RefSeq, Jul 2008]		

— SELECTED CITATIONS —

- **[IF=3.913]** Daijun Wang. et al. MCEMP1 is a potential therapeutic biomarker associated with immune infiltration in advanced gastric cancer microenvironment. GENE. 2022 Oct;840:146760 IHC ;Human. 35905854