

bs-14716R

[Primary Antibody]

CSRNP3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		IHC-P (1:100-500)
GeneID: 80034	SWISS: Q8WYN3	IHC-F (1:100-500)
Target: CSRNP3		IF (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from human FAM130A2: 1-100/585.		ICC/IF (1:100-500)
Purification: affinity purified by Protein A		ELISA (1:5000-10000)
Concentration: 1mg/ml		Reactivity: (predicted: Human, Mouse, Rabbit, Pig, Sheep, Cow, Dog)
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Predicted MW.: 65 kDa
Background: Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.		Subcellular Location: Secreted ,Extracellular matrix