
SLC25A30 Rabbit pAb

Catalog Number: bs-21228R

Target Protein: SLC25A30

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 32 kDa

Entrez Gene: 253512

Swiss Prot: Q5SVS4

Source: KLH conjugated synthetic peptide derived from human SLC25A30: 1-100/291.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Although the outer mitochondrial membrane is permeable to many small metabolites, transport of solutes across the inner mitochondrial membrane is achieved by members of the mitochondrial carrier protein family, such as SLC25A30 (Haguenauer et al., 2005 [PubMed 15809292]).[supplied by OMIM, Mar 2008]